DPV: Dokumentation CGM/FGM

Andreas Hungele

Ramona Ranz

Stand: 20.06.2022 - DPV-Version 7.50.072

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1	Dokumentation in Verlaufsmaske	2
2	Export von CGM-/FGM-Daten aus der jeweiligen Software	4
2.1	Abbott – FreeStyle CoPilot	4
2.2	Abbott – FreeStyle Libre	5
2.3	Abbott - LibreView	5
2.4	Dexcom Clarity	6
2.5	Dexcom Studio	7
2.6	Diasend	7
2.7	Glooko	12
2.8	Medtronic – CareLink Personal	13
2.9	Medtronic – CareLink Pro	13
2.10	myDiabby	14
2.11	Roche - Accu Chek SmartPix Software	15
2.12	Tidepool	15
3	Import von Daten in DPV	17
3.1	Import von CGM-/FGM-Daten	17
3.2	Import von Basalratendaten	19
3.3	Import von Insulindosen	21

Kapitel 1

Dokumentation in Verlaufsmaske

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November 18, 2015

Information update: Dexcom G5™ compatibility with Diasend

We understand that some of the previous information regarding the compatibility between Dexcom G5™ and Diasend has caused some confusion. Therefore, here comes a short clarification, the roll-out plan and an FAQ segment. As mentioned earlier, compatibility will be available in pre-determined markets and the method of data transfer varies.

Data transfer and availability

The new Dexcom G5™ Mobile CGM System has two options. It can be used together with the Dexcom G5™ Mobile Receiver as a “standard CGM unit”. Or it can be used together with an iPhone and the Dexcom app to make the iPhone act as a receiver. This solution is called the Dexcom G5™ Mobile CGM System and the data is wirelessly transferred to the Dexcom cloud.

How to upload data from Dexcom G5™ into diasend®

	US	Rest of the world
Cloud to cloud Dexcom G5™ Mobile CGM System - The data is wirelessly transferred via an iPhone and stored in Dexcom’s cloud. This data is then seamlessly integrated into the diasend® cloud by adding the Dexcom G5™ in the connected apps section in your diasend® Personal account.	Available in Q2 2016*	✓
diasend® Transmitter Data from the Dexcom G5™ Mobile Receiver can be uploaded into diasend® by using a micro USB cable and a diasend® Transmitter at a clinic.	✓	✓
diasend® Uploader Data from the Dexcom G5™ Mobile Receiver can be uploaded into diasend® by using a micro USB cable and the diasend® Uploader software.	Within the next month	Within the next month


* For more information on cloud to cloud availability in the US, please contact Dexcom’s customer service.

How to upload Dexcom G5™ Mobile - cloud to cloud (Currently not available in the US)


Requires an iPhone and Dexcom's G5™ Mobile app

If you are using an iPhone to act as a receiver, it is possible to transfer data wirelessly from Dexcom's cloud directly into diasend®.

- 1 Log into your diasend® Personal account and click on the button **Connect app**.


Connect app

- 2 Click on the **Connect** link for Dexcom G5™ Mobile*.

Connect

* This feature is currently not available in the US. For more information on the cloud to cloud availability in the US, please contact Dexcom's customer service.

- 3 Sign in with your **Dexcom** username and **password** to authorize and add the Dexcom G5™ Mobile to your diasend® Personal account.

FAQ

- 1. When can I upload my Dexcom G5™ Mobile Receiver from home (via my diasend® personal account)?**
 - a. The current time frame is to have an update released within the next month.
 - b. In the meantime, if your clinic has a diasend® Clinic account and a diasend® Transmitter, you can upload at the clinic.
- 2. Will my Dexcom G5™ data sync seamlessly with Diasend?**
 - a. Yes, as long as this is available in your country/region and you have an iPhone. The data will sync to diasend® with a 3 hour delay, due to regulatory restrictions.
- 3. I have an Android phone and/or the Dexcom G5™ Mobile System is not yet available in my country. Can I still upload data via diasend®?**
 - a. Yes. Provided that you have a Dexcom G5™ Mobile Receiver and that your clinic has a diasend® Clinic account and a diasend® Transmitter. If so, you will have the ability to upload at the clinic.
- 4. Can my health care provider upload my Dexcom G5™ data?**
 - a. Yes, as long as your health care provider has a diasend® Clinic account and a diasend® Transmitter, then your device can be uploaded at the clinic.
- 5. Can I view my Dexcom G5™ data from home if my receiver was uploaded at the clinic?**
 - a. Yes. As long as you have entered your Dexcom G5™ serial number into your account, under patient profile.
- 6. My Dexcom G5™ data has wirelessly transmitted to my iPhone, but the data does not show up in my diasend® Personal account, why?**
 - a. There is a 3 hour delay for the data to transfer to diasend due to regulatory restrictions.
- 7. I have an iPhone. How do I seamlessly connect my Dexcom G5™ data to Diasend?**
 - a. Login to your diasend® Personal account. Press the **Connect App** button. Then, under the Dexcom G5™ app, press **connect** and enter your Dexcom credentials.

Any other questions?

Should you have any other questions, please do not hesitate to contact our customer support.

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